



Protocol Booklet

| | |
|------------------------|---|
| Product Code(s) | HB24669 |
| Product Name(s) | BAPTA-AM Janelia Fluor® 549 |
| Purpose | Measurement of intracellular Ca^{2+} in cultured cells |

Please note: This product is for RESEARCH USE ONLY and is not intended for therapeutic or diagnostic use. Not for human or veterinary use



Contents

| | |
|---|---|
| Product Overview | 3 |
| Components & Storage | 3 |
| Protocol | 3 |
| Recipes | 4 |
| BAPTA-AM Janelia Fluor® 549 Loading Solution | 4 |
| HEPES-buffered Hank's Balanced Salt Solution (Assay Buffer) | 4 |
| Guidelines, precautions, troubleshooting | 4 |
| Contact | 4 |
| For customers in the UK, Europe and Rest of the World | 4 |
| For customers in the USA, Canada and South America | 4 |



Product Overview

BAPTA-AM Janelia Fluor® 549 is a membrane permeable, red-shifted (Excitation 546nm, Emission 569nm), intracellular calcium (Ca^{2+}) indicator ($K_d = 310\text{nM}$). It is suitable for measurement of fast calcium dynamics in neurons and cardiomyocytes with excellent photostability and brightness compared to genetically encoded sensors. There are reduced issues with tissue autofluorescence and background fluorescence due to the red-shifted fluorophore. BAPTA-AM Janelia Fluor® 549 is compatible with fluorescence microscopy using TRITC or Cy3 filters.

Components & Storage

BAPTA-AM Janelia Fluor® 549 is provided as either:

| SKU | Component | Quantity | Storage Temperature |
|---------|-----------------------------|---------------------|---------------------|
| HB24669 | BAPTA-AM Janelia Fluor® 549 | 500 μg | -20°C |
| | | 10x50 μg | |

This protocol additionally requires:

| Component | Quantity | Storage Temperature |
|---|------------------|---------------------|
| DMSO | 25 μl | RT |
| Pluronic F-127 | 10mg | 4°C |
| Probenecid | 7.7mg | 4°C |
| Assay Buffer (HEPES-buffered Hank's Balanced Salt Solution (pH = 7.3)*) | 10ml | RT |

* Please see recipe at the end of this protocol book.

Protocol

The following protocol provides general guidelines for using BAPTA-AM Janelia Fluor® 549 to measure intracellular calcium in cultured cells. All loading conditions (dye concentration, temperature, and time) should be optimized for your specific assay, application, and instrumentation.

1. Culture cells following standard protocols to approximately 80-100% confluence.
2. Prepare the loading solution freshly following the below table, vortex well and use within 2 hours.
3. Remove the cell culture medium, briefly wash in plain media (without serum), then add dye loading solution. Recommend volumes are:
 - a. 35mm dish / 6-well plate - 1.5 mL/well,
 - b. 96 well plate - 100 μL /well,
 - c. 384 well plate - 20 μL /well,
4. Incubate in a cell culture incubator at 37°C for 30 minutes.
5. Wash cells twice in fresh assay buffer
6. Read fluorescence using either a plate reader (Excitation: 545nm, Emission 575nm) or image using a fluorescence microscope using a compatible filter set (e.g. Cy3, TRITC).



Recipes

BAPTA-AM Janelia Fluor® 549 Loading Solution

| Component | Concentration | Quantity | Notes |
|-----------------------------|---------------|----------|--|
| BAPTA-AM Janelia Fluor® 549 | 4.1µM | 50µg | Dissolve in DMSO then aliquot and store any unused dye at -20°C |
| Assay Buffer | 1X | 10ml | Normally HEPES buffered HBSS but other buffers have been also successfully used. |
| Pluronic F-127 | 0.1% | 10mg | Surfactant that helps the dissolution of dye therefore ensuring even dye distribution and cellular loading. |
| Probenecid | 2.7mM | 7.7mg | Anion transport inhibitor that improves intracellular dye retention. Not required for all cell types, it is recommended in most cases to optimize assay performance. |

Please note: Combine components then vortex thoroughly. Use within 2 hours of creation. Do not freeze.

HEPES-buffered Hank's Balanced Salt Solution (Assay Buffer)

| Component | MW (g/mol) | g/L | Concentration (mM) |
|--------------------------------|------------|-------|--------------------|
| Calcium Chloride | 110.98 | 0.14 | 1.26 |
| Magnesium Chloride Hexahydrate | 203.30 | 0.1 | 0.49 |
| Magnesium Sulfate Heptahydrate | 246.47 | 0.1 | 0.41 |
| Potassium Chloride | 74.55 | 0.4 | 5.33 |
| Potassium Phosphate Monobasic | 136.09 | 0.06 | 0.44 |
| Sodium Bicarbonate | 84.01 | 0.35 | 4.17 |
| Sodium Chloride | 58.44 | 8 | 138.00 |
| Sodium Phosphate Dibasic | 141.96 | 0.048 | 0.34 |
| D-Glucose (Dextrose) | 180.16 | 1 | 5.56 |
| HEPES | 238.30 | 4.76 | 20.00 |

Please note: Add all components to dH₂O, mix well then adjust to pH 7.3

Guidelines, precautions, troubleshooting

Please contact our technical support team at technicalhelp@hellobio.com for advice on how to resolve any problems encountered when using this product. Observe safe laboratory practice and consult the safety datasheet. Please see the datasheet on our website for general guidelines, precautions, limitations on the use of the product.

Contact

For customers in the UK, Europe and Rest of the World

Customer Care customercare@hellobio.com

Technical support technicalhelp@hellobio.com

By telephone: +44(0)117 318 0505

By fax: +44(0)117 981 1601

For customers in the USA, Canada and South America

Customer Care customercare-usa@hellobio.com

Technical support technicalhelp@hellobio.com

By telephone: +1-609-683-7500

By fax: +1-609-228-4994



Opening hours: 8.30 am - 5.00 pm GMT weekdays

Opening
hours:

9.00 am - 5.00 pm EST weekdays