



Frequently Asked Questions

- **What does the product look like?**

- Appearance:

- The products are supplied as a dried film of blue material adhering to the inside of the conical bottom of a 'champagne flute' brown glass vial
 - The product should be visible as a dark smear through the side of the vial


- **The fluorophore is red, but the solution is blue - have I made it up properly?**

- Reconstitution:

- If dissolved in DMSO at 10mM according to the instructions in the datasheet, the product will appear to be a dark blue solution – this is correct
 - For sub-aliquoting in slightly larger volumes, the product can also be reconstituted at 1mM by adding 10x the volume of DMSO stated
 - The product can be further diluted in DMSO to eg 100µM and will retain its blue colour
 - Final dilution to ~ 100nM is in buffer (+ 0.1% DMSO)

- **What excitation and emission wavelengths and filter sets should be used?**

- The BODIPY-630/650 fluorophore has its excitation peak at 625nm and its emission peak at 641nm
- These peaks are close, so narrow band-pass filters are needed
- For excitation, use a narrow-bandpass filter with a peak at 630nm
 - For laser-based systems the 633nm ruby laser is ideal for excitation
- For emission, use a narrow-band filter with a peak at 650nm
- Shorter wavelength excitation and longer wavelength emission filters can be used but signals will be weaker, and there may be cross-talk with other fluorophores



- **Cells don't bind the fluorescent ligand**

- Likely causes:

- Cells are not expressing the receptor

- Compare with cells from another source
 - Compare with cells at different growth stages
 - Use the ligand for FACS to select cells expressing the receptor

- The receptor is not reaching the cell membrane

- Ligands do not readily cross the cell membrane, so this is difficult to test
 - An antibody approach with fixed, permeabilised cells might confirm this


- The receptor is unable to bind the ligand

- Fixing the cells can prevent ligand binding – use live cells
 - Too much solvent (DMSO) in buffers can damage the cells

- **Cells label everywhere (cytoplasm, nucleus), not just the membrane**

- This is likely to be because the cells are not healthy and have taken up the ligand non-specifically.

- Check the culture conditions – cells should not be over-confluent
- Check the assay conditions – pH should be 7.4 in isotonic buffer (eg. HBS or PBS)
- Check there is minimal solvent used (eg DMSO < 0.1%)
- Check the cells have not been fixed and permeabilised



- Even in a healthy culture, there will be dead cells that show non-specific staining.


- Make sure cell-plates are handled carefully and removed from the incubator immediately before use, not left sitting on the bench for long periods

- **Control (non target-expressing) cells label as well as the target cells**

- Likely cause: high non-specific binding to control and target cells
 - Ensure both control and target cell lines are healthy (see previous FAQ)
 - Ensure the ligand is not being used at too high concentration (maximum should be 100nM)
- Endogenous receptor expression by the host cell line
 - Check if the ligand is binding to an endogenous receptor expressed by the host cells by blocking with an appropriate unlabelled competitor
- Cross-reaction to other receptors
 - the ligands may bind to other receptors
 - Use the pharmacology of the unlabelled drug from which the ligand is derived as a guide to potential additional targets to which the ligand may be binding
 - Test this with selective unlabelled blockers for the off-target receptor

- **There is very high background fluorescence**

- BODIPY-630/650 should not show this because it's fluorescence is quenched in aqueous solution
- If this is reported, it is likely to be due to a high concentration of the ligand in solution around the cells
 - Using the ligand at a lower final concentration (30 – 100nM gives the best signal : background).
 - Add a brief (10 min), gentle wash step to reduce the background fluorescence



- **Cells label well, and show a good saturation-binding curve**

- but the maximal saturation binding level is lower than expected

- This may be due to the fluorescence detector reaching saturation, rather than the cells reaching binding saturation.

- Try reducing the gain on the detector
 - Try taking measurements on an alternative instrument

